EL30

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Revisão 00 - Jul/2026

NKL®

www.nkl.com.br

Tabela de Memória: T26

A — Tempo Total: 30min

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FARMER, A. D.; STRZELCZYK, A.; FINISGUERRA, A.; GOURINE, A. V.; GHARABAGHI, A.; HASAN, A.; et al. International consensus based review and recommendations for minimum reporting standards in research on transcutaneous vagus nerve stimulation (Version 2020). Frontiers in Human Neuroscience, v. 14, p. 568051, 2021.

YAP, J. Y. Y.; KEATCH, C.; LAMBERT, E.; WOODS, W.; STODDART, P. R.; KAMENEVA, T. Critical review of transcutaneous vagus nerve stimulation: challenges for translation to clinical practice. Frontiers in Neuroscience, v. 14, p. 284, 2020.

BONAZ, B.; SINNIGER, V.; PELLISSIER, S. Vagus nerve stimulation: a new promising therapeutic tool in inflammatory bowel disease. Journal of Internal Medicine, v. 282, n. 1, p. 46–63, 2017.

TRACEY, K. J. The inflammatory reflex. Nature, v. 420, p. 853–859, 2002.

TRACEY, K. J. Reflex control of immunity. Nature Reviews Immunology, v. 9, p. 418–428, 2009.

BREIT, S.; KUPFERBERG, A.; ROGLER, G.; HASLER, G. Vagus nerve as modulator of the brain–gut axis in psychiatric and inflammatory disorders. Frontiers in Psychiatry, v. 9, p. 44, 2018.

FRANGOS, E.; ELLRICH, J.; KOMISARUK, B. R. Non-invasive access to the vagus nerve through the external ear: fMRI evidence in humans. Brain Stimulation, v. 8, n. 3, p. 624–636, 2015.

YAKUNINA, N.; KIM, S. S.; NAM, E. C. Optimization of transcutaneous vagus nerve stimulation using functional MRI. Neuromodulation, v. 20, n. 3, p. 290–300, 2017.

REDGRAVE, J. N.; DAY, D.; LEUNG, H.; LAUD, P. J.; ALI, A.; LINDERT, R.; et al. Safety and tolerability of transcutaneous vagus nerve stimulation in humans: a systematic review. Brain Stimulation, v. 11, n. 6, p. 1225–1238, 2018.

GEORGE, M. S.; ASTON-JONES, G. Noninvasive techniques for probing neurocircuitry: vagus nerve stimulation. Biological Psychiatry, 2010.

GROVES, D. A.; BROWN, V. J. Vagal nerve stimulation: a review of its applications and potential mechanisms. Neuroscience & Biobehavioral Reviews, v. 29, p. 493–500, 2005.

ENGINEER, N. D.; KIMBERLEY, T. J.; PRUDENTE, C. N.; DAWSON, J.; TARVER, W. B.; HAYS, S. A. Targeted vagus nerve stimulation for rehabilitation after stroke. Frontiers in Neuroscience, v. 13, p. 280, 2019.

DAWSON, J.; LIU, C. Y.; FRANCISCO, G. E.; CRAMER, S. C.; WOLF, S. L.; DIXIT, A.; et al. Vagus nerve stimulation paired with rehabilitation for upper limb motor function after ischemic stroke. The Lancet, 2021.

KIMBERLEY, T. J.; PIERCE, D.; PRUDENTE, C. N.; FRANCISCO, G. E.; YOZBATIRAN, N.; SMITH, P.; et al. Vagus nerve stimulation paired with upper limb rehabilitation after chronic stroke. Stroke, 2018.

CLARK, G. T. Clinical electrophysiology: electrotherapy and electrophysiologic testing. 4. ed. Wolters Kluwer, 2021.

ROBINSON, A. J.; SNYDER-MACKLER, L. Clinical electrophysiology: electrotherapy and electrophysiologic testing. 4. ed. Lippincott Williams & Wilkins, 2008.

CAMERON, M. H. Physical agents in rehabilitation. 5. ed. Elsevier, 2017.

KITCHEN, S.; BAZIN, S. Clayton's electrotherapy. 12. ed. Elsevier, 2020.

WATSON, T. Electrotherapy: evidence-based practice. Elsevier, 2008.

SHEFFLER, L. R.; CHAE, J. Neuromuscular electrical stimulation in neurorehabilitation. Muscle & Nerve, v. 35, p. 562–590, 2007.

POPOVIC, D. B.; SINKJAER, T. Control of movement for the physically disabled. Springer.

CARNABY, G. D.; CRARY, M. A. Examining the evidence on neuromuscular electrical stimulation for swallowing. Perspectives of the ASHA Special Interest Groups.

FREED, M. L.; FREED, L.; CHATBURN, R. L.; CHRISTIAN, M. Electrical stimulation for swallowing disorders caused by stroke. Respiratory Care, v. 46, p. 466–474, 2001.

LUDLOW, C. L.; HUMBERT, I.; SAXON, K.; POLETTO, C.; SONIES, B.; CRUJIDO, L. Effects of surface electrical stimulation both at rest and during swallowing in chronic pharyngeal dysphagia. Dysphagia, v. 22, p. 1–10, 2007.

SPEYER, R.; BAIJENS, L.; HEIJNEN, M.; ZWIJNENBERG, I. Effects of therapy in oropharyngeal dysphagia by means of neuromuscular electrical stimulation: a systematic review. Dysphagia, 2010.

BATH, P. M.; LEE, H. S.; EVERTON, L. F. Swallowing therapy for dysphagia in acute and subacute stroke. Cochrane Database of Systematic Reviews, 2018.

RESTIVO, D. A.; HAMDY, S. Pharyngeal electrical stimulation for dysphagia. Neurotherapeutics, 2018.

HAMDY, S.; ROTHWELL, J. C.; AZIZ, Q.; THOMPSON, D. G. Long-term reorganization of human motor cortex driven by short-term sensory stimulation. Nature Neuroscience, 1998.

HUCKABEE, M. L.; STEELE, C. M. An analysis of surface electromyography in swallowing research. Journal of Speech, Language, and Hearing Research.

BEHRMAN, A. Speech and voice science. Plural Publishing.

BOONE, D. R.; McFARLANE, S. C.; VON BERG, S. L.; ZRAICK, R. I. The voice and voice therapy. 10. ed.

STEMPLE, J. C.; ROY, N.; KLABEN, B. K. Clinical voice pathology. 6. ed.

DWORKIN, J. P.; MELECA, R. J. Laryngeal electromyography.

ASHA – AMERICAN SPEECH-LANGUAGE-HEARING ASSOCIATION. Adult dysphagia practice portal.

EUROPEAN SOCIETY FOR SWALLOWING DISORDERS (ESSD). Clinical practice guidelines for oropharyngeal dysphagia.

CÔTÉ, M. P.; MURRAY, L. M.; KNIKOU, M. Spinal control of locomotion: individual neurons, their circuits and functions. Nature Reviews Neuroscience.

CRUCCU, G.; AZIZ, T. Z.; GARCIA-LARREA, L.; HANSSON, P.; JENSEN, T. S.; LEFAUCHEUR, J. P.; et al. EFNS guidelines on neurostimulation therapy for neuropathic pain. European Journal of Neurology, 2007.

LEFAUCHEUR, J. P.; ALEMAN, A.; BAEKEN, C.; BENNINGER, D. H.; BRUNELIN, J.; DI LAZZARO, V.; et al. Evidence-based guidelines on the therapeutic use of transcranial electrical stimulation. Clinical Neurophysiology, 2017.

WORLD HEALTH ORGANIZATION. WHO guidelines on rehabilitation. Geneva.

MELZACK, R.; WALL, P. D. Pain mechanisms: a new theory. Science, v. 150, n. 3699, p. 971–979, 1965.

JOHNSON, M. I. Transcutaneous electrical nerve stimulation (TENS): research to support clinical practice. Oxford University Press, 2014.

SLUKA, K. A.; WALSH, D. Transcutaneous electrical nerve stimulation: basic science mechanisms and clinical effectiveness. Journal of Pain, v. 4, n. 3, p. 109–121, 2003.

VANCE, C. G. T.; DAILEY, D. L.; RAKEL, B. A.; SLUKA, K. A. Using TENS for pain control: the state of the evidence. Pain Management, v. 4, n. 3, p. 197–209, 2014.

CRUCCU, G.; TRUINI, A. Refractory trigeminal neuralgia: non-surgical treatment options. CNS Drugs, v. 27, n. 2, p. 91–96, 2013.

SCHIFFMAN, E.; OHRBACH, R.; TRUELOVE, E.; et al. Diagnostic criteria for temporomandibular disorders (DC/TMD). Journal of Oral & Facial Pain and Headache, v. 28, n. 1, p. 6–27, 2014.

ROBINSON, A. J.; SNYDER-MACKLER, L. Clinical electrophysiology: electrotherapy and electrophysiologic testing. 4. ed. Philadelphia: Wolters Kluwer/Lippincott Williams & Wilkins, 2021.

WATSON, T.; NUSSBAUM, E. L. Electrophysical modalities: evidence-based practice. 13. ed. Elsevier, 2020.

CAMERON, M. H. Physical agents in rehabilitation: an evidence-based approach to practice. 6. ed. Elsevier, 2023.

KITCHEN, S.; BAZIN, S.; CLAYTON, K. Clayton's electrotherapy: theory and practice. 12. ed. Elsevier, 2022.

JOHNSON, M. I. Transcutaneous electrical nerve stimulation (TENS): research to support clinical practice. Oxford University Press, 2014.

SLUKA, K. A.; BJORDAL, J. M.; MARCHAND, S.; RAKEL, B. A. What makes transcutaneous electrical nerve stimulation work? Making sense of the mixed results in the clinical literature. Physical Therapy, v. 93, n. 10, p. 1397–1402, 2013.

VANCE, C. G. T.; DAILEY, D. L.; RAKEL, B. A.; SLUKA, K. A. Using TENS for pain control: the state of the evidence. Pain Management, v. 4, n. 3, p. 197–209, 2014.

JOHNSON, M. I.; PALEY, C. A.; JONES, G.; et al. Efficacy and safety of transcutaneous electrical nerve stimulation (TENS) for acute and chronic pain: overview of systematic reviews. BMJ Open, 2022.

JOHNSTON, T. E.; KELLER, S.; DENZER-WEILER, C.; BROWN, L. A clinical practice guideline for the use of ankle-foot orthoses and functional electrical stimulation post-stroke. Journal of Neurologic Physical Therapy, v. 45, p. 112–196, 2021.

ASSOCIATION OF CHARTERED PHYSIOTHERAPISTS IN NEUROLOGY (ACPIN). Evidence based clinical practice guidelines for the use of functional electrical stimulation to improve mobility in adults with lower limb impairment due to an upper motor neuron lesion. 2022.

SHEFFLER, L. R.; CHAE, J. Neuromuscular electrical stimulation in neurorehabilitation. Muscle & Nerve, v. 35, n. 5, p. 562–590, 2007.

POPOVIĆ, D. B. Advances in functional electrical stimulation (FES). Journal of Electromyography and Kinesiology, v. 24, n. 6, p. 795–802, 2014.

HOUGHTON, P. E.; NUSSBAUM, E. L.; HOENS, A. M.; et al. Electrophysical agents: contraindications and precautions. Physiotherapy Canada, 2010.

SILVA, F. P.; et al. Electrophysical agents in clinical practice of orthopedic and sports physiotherapy. Fisioterapia em Movimento, 2020.